

# High efficiency in bacteria retention

Nylon 66 membrane

Pre-flushed with DI non pyrogenic water


Individually integrity tested before packing

## ABSO NY

### N66 nylon for enhanced bio-reduction performance and submicron particle removal

ABSO NY cartridges are made of inherently hydrophilic N66 membrane, specifically designed for bio-reduction and final filtration of pharma solutions. The N66 membrane has an uniform pore distribution, repeatedly from batch to batch and the efficiency can be verified through integrity test.

The use of thermo-bonding and ultrasound bonding processes in assembly without the use of resins or adhesives reduces the level of extractables and makes the whole component compact and resilient. Manufactured in clean room and individually tested to verify integrity before packing.



## > ABSO NY

100% inherently hydrophilic,  
no surfactants or wetting  
agents

## TECHNICAL SPECIFICATIONS

- N66 membrane, uniform filtration across the entire surface, non-fiber migration
- Absolute retention ratings, high margins of operational safety
- pre-flushed with DI non pyrogenic water (<0.25 EU/ml)
- Manufactured in clean room
- All materials meet FDA CFR 21; Plastics compliant with USP Class VI “Plastic Biosafety”
- In compliance with the EU regulation for food contact in EC Countries (10/2011 + amendments; 1935/2004;1895/2005)
- 100% tested for integrity prior to packing

## OPERATING CONDITIONS

<b>Max operating pressure (<math>\Delta p</math>)</b>	80°C @ 1.0 bar 20°C @ 5.0 bar
<b>Recommended replacement pressure drop</b>	2.0 bar @20°C
<b>Suggested operating pressure range</b>	0.1 to 1.0 bar

## PORE SIZE RATING & TYPICAL DATA

0.2um – 0.45um – 0.65um

<b>Membrane Retention Efficiency - 350l/h/10”</b>	<b>Microbial Titer Reduction (<math>T_R</math>) ASTM F 838-05</b>
>99.99 %	$T_R > 10^7$ cm <sup>2</sup> // 0.2um [Brevundimonas Diminuta]
>99.99 %	$T_R > 10^7$ /cm <sup>2</sup> // 0.45um – 0.65um [Serratia Marcescens]

## FLOW RATE

<b>Water Flow 20°C @0.1 bar /10”</b>	0.20um	0.45 um	0.65um
<b>Typical Flow Rate</b>	7 L/min	12 L/min	15L/min

Extrapolation for multiple housings and higher pressure drops is acceptable, but as flows increase the pressure drop of the housing becomes more apparent

## FILTRATION AREA

0.6m<sup>2</sup>/10”

## SANITIZATION

Steam : 125°C/30mins – Hot Water : 80°C/30mins – Autoclave : 121°C/30mins

## MATERIALS OF CONSTRUCTION

<b>Filtering media</b>	Nylon 66
<b>Supports</b>	Polypropylene
<b>Inner sleeve</b>	Polypropylene + ss316 insert
<b>Connections and tip</b>	Polypropylene
<b>Gaskets</b>	Silicone (standard), EPDM, Viton, FEP

## SEALING

Ultrasound / Heat sealing

## DIMENSION

<b>Length</b>	254mm (10”)–508mm (20”)–762mm (30”)–1016mm (40”)
<b>Outer diameter</b>	69mm
<b>Inner diameter</b>	26mm

## TRACEABILITY

Each filter element is identified by a lot number for complete traceability.